

VACCI^XCELL

ADHERENT BIOPROCESSING SPECIALIST



West Coast
Representative

Margaret Stava
562 355-5859
margaret@biopharmequip.com

Cell Processing Isolator

Highly Adaptable Solution for Aseptic Bioprocessing

The Esco VacciXcell Cell Processing Isolator is a one-stop solution which integrates various cell processing equipment in a cGMP isolator necessary for aseptic bioprocessing.



Key Features:

- ✓ Designed in accordance with international cGMP standards
- ✓ Full unidirectional airflow with H14 filters providing ISO Class 5 environment
- ✓ Class 2 containment enclosure as per ISO10648-2 standards
- ✓ Matched environment with the incubation chamber
- ✓ PLC-based monitoring and control system

Equipment Integration:

- ✓ CO₂ Incubator
- ✓ Centrifuge
- ✓ Microscope
- ✓ BioVap™ decontamination system
- ✓ Glove Leak Tester
- ✓ Particle Monitoring