




Safeport TransferBag (Patent pending)

With the SafePort Bag, we offer a cost-effective oneway option for a secure transfer with OEB 5 (<math><1 \mu\text{g per m}^3</math>) containment-level.

With the development of the SafePort Bag, Lugaia Containment Solutions has placed the highest priority on easy handling, high working safety, universal usability and quality.

The design of the SafePort Bag avoids product loss, thus ensuring a high yield of your products.

- Clean room manufacturing according to the state of the art guidelines with USP / EP certificates.
- Standard version manufactured in the cleanroom classes ISO 8 and ISO 7.
- Special version with sampling sleeve and customized films available.



Standard sizes

Volume (l)	Bag width (mm)	Bag length (mm)
5	300	1100
25	500	1300
40	600	1500

Special sizes on request!

Connection options:

- Inlet and Outlet SafePort
- Inlet O-Ring, Outlet SafePort
- Inlet open, Outlet SafePort
- InletTriClamp, Outlet SafePort

