





Overview of standard types

The Lugaia TransferBags meet the high requirements of the pharmaceutical and biotech sector. Standard versions are available in a permanently antistatic film. If desired, the bags can be made from a customer-specific film in the clean room. The design of the TriClamps guarantees a complete emptying of the bag and thus ensures the maximum product yield.



| Volume (l) | Bag-width W1 (mm) | Bag-length L (mm) | TriClamp-connection (*) | Ø-Connection for O-Ring (mm) |
|------------|-------------------|-------------------|-------------------------|------------------------------|
| 5 | 400 | 400 | 4 | 170 |
| 10 | 400 | 700 | 4/6 | 170 |
| 15 | 400 | 1100 | 4/6 | 170 |
| 20 | 500 | 1100 | 4/6 | 230 |
| 25 | 500 | 1400 | 4/6 | 230 |
| 30 | 700 | 1200 | 4/6 | 230 |
| 40 | 700 | 1600 | 4/6 | 230 |

Contact:
info@biopharmequip.com
 562 355-5859

www.lugaia.com

Lugaia AG Containment Solutions Industriezone Basper 31, 3942 Raron, Schweiz
 Telefon: +41 (0) 27 948 45 55, Telefax: +41 (0) 27 948 45 50, anfrage@lugaia.ch

Lugaia Deutschland GmbH Marie-Curie-Straße 6, 85055 Ingolstadt, Deutschland
 Telefon: +49 (0) 841 31 96 17 90, Telefax: +49 (0) 841 31 96 17 95, info@lugaia.de